

# LapMan<sup>®</sup> Dynamic Laparoscope Driver



*Get a stable image*

*Control your vision yourself*

*Make better use of your staff*

**“Control your vision yourself”**

**Med**  
SYS

[www.medsys.be](http://www.medsys.be)

## Introducing the LapMan<sup>®</sup> into the OR will :

Give back the surgeon the control of his vision

Improve image stability

Produce less eye fatigue

Free the Assistant for more rewarding tasks

Allow Solo Surgery

### Features :

Rolling Base makes positioning easy, on either side of the operating table

Connection to laparoscope is easy, fast and scope independent

Set up takes less than 5 minutes

### LapStick<sup>®</sup> :

Fits to any laparoscopic instrument (tailor-made clips)

Surgeon's index controls all scope movements

Sterile Disposable Device



#### Technical specifications

##### LapMan<sup>®</sup>

Weight : 50 Kgs

Dimensions :

Height : 1,2 - 1,5 m

Length : 0,7 m

Width : 0,36 m

##### LapStick<sup>®</sup>

Weight : 35 gr

Dimensions :

8,5 cm x 3 cm x 2 cm

**“Control your vision yourself”**

Medsys s.a. • 39, rue Chainisse  
B-5030 Gembloux • Belgium  
Tel. : +32 81 56 84 84 • Fax : +32 81 56 88 84  
E-mail : [info@medsys.be](mailto:info@medsys.be)  
<http://www.medsys.be>

**Med**  
**SYS**